

VETERINARY INSTRUMENTATION

Broadfield Road, Sheffield UK. S8 0XL
 Tel: 0845 130 9596
 Fax: 0845 130 8687
 info@vetinst.com www.vetinst.com

CCL SELECTION OF NYLON

WEIGHT KGS (LBS)	STRENGTH OF NYLON	NYLON	NYLON + CRIMP	NYLON, CRIMP + SWAGED ON FABELLA NEEDLE	NYLON + 1 x SWAGED ON FABELLA NEEDLE - DOUBLE LENGTH	NYLON + 2 x CRIMP + 1 x SWAGED ON FABELLA NEEDLE - DOUBLE LENGTH
UP TO 15 (33)	50 lbs	091141 (500mm)	091142 (500mm)	091155 (500mm)	091160 (1000mm loop)	091165 (1000mm loop)
20-25 (44-55)	80 lbs	091139 (800mm)	091137 (800mm)	091156 (800mm)	091161 (1600mm loop)	091166 (1600mm loop)
30 (66)	100 lbs	091148 (800mm)	091147 (800mm)	091157 (800mm)	091162 (1600mm loop)	091167 (1600mm loop)
VERY LARGE DOGS	2 X 100lbs	2 x 091148	2 x 091147	2 x 091157	1 x 091162	1 x 091167

All packs above are supplied pre-sterilised.

Suggested nylon usage is for guidance only. The surgeon must make their decision based on the clinical case taking into account obesity, activity levels etc. If in doubt use the maximum. The nylon will always break after around 6 weeks if it is providing support. The break is usually at the fabella which appears to be the point of maximum abrasion. Breakage at this time does not affect the outcome.