











# Dental Burs - Many New Shapes

To aid correct bur identification, use the following steps.

1. Identify the head shape (Green Column)
2. Select the diameter of head (White Column).
3. Select the shaft type (Top of table, Pink, Yellow, Blue, Purple)
4. Select Carbide or Diamond

Head Shape	Diameter	FG 		RA 		HP 			
		Carbide 5 Pack	Diamond 5 Pack	Carbide 5 Pack	Diamond 5 Pack	Carbide Singles	Carbide 5 Pack	Diamond Singles	Diamond 5 Pack
Small Round 	0.5	5FG 1/4	5FG 1/4D						
Medium Round	1.4	5FG4	5FG4D	5RA4	5RA4D	HP4	5HP4	HP4D	5HP4D
Large Round	1.8	5FG6	5FG6D	5RA6	5RA6D	HP6	5HP6	HP6D	5HP6D
Extra Large Round	2.1	5FG8	5FG8D	5RA8	5RA8D	HP8*	5HP8	HP8D	5HP8D
Extra Small Taper Fissure 	0.8	5FG698		5RA698		HP698	5HP698		
Small Taper Fissure	1.2	5FG701		5RA170L	5RA170LD	HP701	5HP701		
Medium Taper Fissure	1.6	5FG702	5FG702D			HP702	5HP702		
Medium Taper Long Head	1.6	5FG702L	5FG702LD	5RA702L	5RA702LD	HP702L	5HP702L		
Large Taper Fissure	2.1					HP703	5HP703	HP703D	5HP703D
Extra Large Taper Fissure	2.3					HP704*	5HP704	HP704D	5HP704D
Extra Small Flat Fissure 	0.9	5FG556		5RA556		HP556	5HP556		
Small Flat Fissure	1	5FG557		5RA557L	5RA557LD	HP557L	5HP557L	HP557D	5HP557LD
Medium Flat Fissure	1.6	5FG5660		5RA560		HP560	5HP560		
Large Flat Fissure	2.3					HP563	5HP563	HP563D	5HP563D
Small Inverted Cone 	0.6	5FG331/2							
Medium Inverted Cone	1	5FG35	5FG35D						
Large Inverted Cone	1.4	5FG37		5RA37	5RA37D				
Extra Large Inverted Cone	1.8					HP39	5HP39	HP39D	5HP39D
Extra Small Pear 	0.6	5FG329							
Small Pear	0.8	5FG3330							
Medium Pear	1	5FG331							
Large Pear	1.2	5FG332				HP332	5HP332		
Medium Domed Fissure 	1	5FG1557							
Large Domed Fissure	1.2	5FG1558	5FG1558D						

Additional Burs Extra Long	Diameter	HL 			
		Carbide Singles	Carbide 5 Pack	Diamond Singles	Diamond 5 Pack
Extra Large Taper Fissure		HL704	5HL704	HL704D	5HL704D
Large Flat Fissure		HL563*	5HL563	HL537D	5HL537D

DENTAL BUR SINGLE PRICES		
BUR PREFIX	CARBIDE	DIAMOND
HP	£3.50	£4.50
HL	£4.50	£5.25

DENTAL BUR PRICES (5 PACK)		
BUR PREFIX	CARBIDE	DIAMOND
5FG	£13.75	£16.50
5RA	£13.75	£16.50
5HP	£13.75	£16.50
5HL	£19.75	£23.75

\* in Rabbit bur Kit