

VETERINARY INSTRUMENTATION

Broadfield Road, Sheffield UK. S8 0XL

Tel: 0845 130 9596

Fax: 0845 130 8687

info@vetinst.com www.vetinst.com

LIGAFIBA QUICK REFERENCE CHART

LF150500	LIGAFIBA 150LB X 500MM (STERILE)	<15KG
LF250600	LIGAFIBA 250LB X 600MM (STERILE)	<50KG
LF500600	LIGAFIBA 500LB X 600MM (STERILE)	>50KG
LFCRIMP150	CRIMP TUBE FOR 150LB LIGAFIBA (NON STERILE)	N/A
LFCRIMP250	CRIMP TUBE FOR 250LB LIGAFIBA (NON STERILE)	N/A
LFCRIMP500	CRIMP TUBE FOR 500LB LIGAFIBA (NON STERILE)	N/A
HIPTOGLF250	250LB LIGAFIBA HIP TOGGLE (250LB)	<35KG
LFIS253518SS	LIGAFIBA INTERFERENCE SCREW 3.5 DIA 18MM LONG STAINLESS	250LB LIGAFIBA
LFIS253518TI	LIGAFIBA INTERFERENCE SCREW 3.5MM DIA 18MM LONG TITANIUM	250LB LIGAFIBA
LFIS354523SS	LIGAFIBA INTERFERENCE SCREW 4.5 DIA 23MM LONG STAINLESS	500LB LIGAFIBA
LFIS354523TI	LIGAFIBA INTERFERENCE SCREW 4.5MM DIA 23MM LONG TITANIUM	500LB LIGAFIBA
LFISOKIT	LIGAFIBA ISO TOGGLE STARTER KIT (250LB + 500LB)	
LFITS150	LIGAFIBA ISO TOGGLE SUTURE 150LB USES 2MM STD DRILL	15KG
LFITS250	LIGAFIBA ISO TOGGLE SUTURE 250LB USES 2.5MM STD DRILL	15-40KG
LFITS500	LIGAFIBA ISO TOGGLE SUTURE 500LB USES 3.5MM STD DRILL	>40KG
LFLS150	LIGAFIBA LATERAL SUTURE 150LB (STR & CVD NEEDLE)	<15KG
LFLS250	LIGAFIBA LATERAL SUTURE 250LB (STR & CVD NEEDLE)	<50KG
LFLS500	LIGAFIBA LATERAL SUTURE 500LB (STR & CVD NEEDLE)	>50KG
LFIS140TC	LIGAFIBA SCISSOR 140MM 5.5" COF T.C.	NA
	LIGAFIBA TENDON REPAIR SUTURE	NA

