

Pin Reference Chart

		CLAMP STYLE		BAR DIAMETER		Small		Medium		Large		
		KE Plus	Small	1/8" / 3.2mm	Medium 3/16" / 4.8mm	Large 5/16" / 8mm	Standard	Mini 1/8" / 3.2mm	Small 1/4" / 6.3mm	Large 1/3" / 9.5mm	Linear	Medium
TROCHAR PINS												
001556	Trochar Pin 2mm	●		●				●				
001557	Trochar Pin 3mm	●		●		●		●		●		
001558	Trochar Pin 4mm			●		●				●	●	●
CORTICAL NEGATIVE THREADED PINS (ELLIS TYPE)												
001528	Ellis pin 1.1m	●						●				
001529	Ellis pin 1.6mm	●		●				●				
001550	Ellis pins 2mm	●		●		●		●			●	
001530	Ellis pins 3/32" - (2.4mm)	●		●		●		●			●	
001550A	Ellis pins 7/64" - (2.7mm)	●		●		●		●			●	
001551	Ellis pins 3mm	●		●		●		●		●	●	
001531	Ellis pins 1/8" - (3.2mm)	●		●		●		●		●	●	
001551A	Ellis pins 9/64" - (3.6mm)	●		●		●		●		●	●	●
001552	Ellis pins 4mm							●		●	●	●
BICORTICAL FINE NEGATIVE THREADED PINS												
Code	Shaft Diameter											
BCNET08	0.8mm									●		
BCNET10	1.0mm	●						●				
BCNET11	1.1mm	●		●				●				
BCNET15	1.5mm	●		●				●				
BCNET16	1.6mm	●		●		●		●			●	
BCNET20	2.0mm	●		●		●		●			●	
BCNET24	2.4mm	●		●		●		●		●	●	
BCNET30	3.0mm	●		●		●		●		●	●	●
BCNET35	3.5mm	●		●		●		●		●	●	●
BCNET40	4.0mm			●		●				●	●	●
BCNET50	5.0mm			●		●					●	●
POSITIVE THREADED PINS - END THREAD												
CET0018	1.6mm shank 2.0mm threads	●		●				●			●	
CET0019	2.0mm shank 2.4mm threads	●		●		●		●			●	
CET0020	2.4mm shank 3.2mm threads	●		●		●		●			●	
CET0020A	2.7mm shank 3.5mm threads	●		●		●		●			●	
CET0020B	3.0mm shank 3.6mm threads	●		●		●		●		●	●	
CET0021	3.2mm shank 4.0mm threads	●		●		●		●		●	●	
CET0024A	3.5mm shank 4.3mm threads			●		●		●		●	●	●
CET0022	4.0mm shank 4.8mm threads			●		●		●		●	●	●
POSITIVE THREADED PINS - MID THREAD												
CMT0018	1.6mm shank 2.0mm threads	●		●				●			●	
CMT0019	2.0mm shank 2.4mm threads	●		●		●		●			●	
CMT0023	2.4mm shank 3.2mm threads	●		●		●		●			●	
CMT0023A	2.7mm shank 3.5mm threads	●		●		●		●			●	
CMT0023B	3.0mm shank 3.6mm threads	●		●		●		●		●	●	
CMT0024	3.2mm shank 4.0mm threads	●		●		●		●		●	●	●
CMT0024A	3.5mm shank 4.3mm threads			●		●		●		●	●	●
CMT0025	4.0mm shank 4.8mm threads			●		●		●		●	●	●
POSITIVE THREADED PINS - CANCELLOUS												
CANET0026	2.4mm shank 3.5mm threads	●		●		●		●			●	
CANET0027	3.2mm shank 4.8mm threads	●		●		●		●		●	●	
CANMT0029	2.4mm shank 3.5mm threads	●		●		●		●		●	●	
CANMT0030	3.2mm shank 4.8mm threads			●		●				●	●	

VETERINARY INSTRUMENTATION

Broadfield Road, Sheffield UK. S8 0XL

Tel: 0845 130 9596

Fax: 0845 130 8687

info@vetinst.com www.vetinst.com

Positive Pins Insertion Data

Pin Codes	Shank Size (MM)	Thread Size (MM)	KE+ Clamp	Drill bit code	Pre drill size
CET/ CMT0018	1.6	2.0	S	H090101	1.5
CET/ CMT0019	2.0	2.4	S	H090202	1.9
CET0020/ CMT0023	2.4	3.2	S/M	H090203	2.3
CET0020A/ CMT0023A	2.7	3.5	M	H090112	2.5
CET0020B/ CMT0023B	3.0	3.6	M/L	H090104	2.7
CET0021/ CMT0024	3.2	4.0	M/L	H090205	3.0
CET/ CMT0024A	3.5	4.3	M/L	H090105S H090105L	3.2
CET0022/ CMT0025	4	4.8	M/L	H090207	3.8

Negative Pins (Ellis Type and Bicortical -VE Threaded) Insertion Data

Pin Code Ellis Bicort	Shank Size (MM)	KE+ Clamp	Pre Drill Size	Drill bit code
001528 BCNET11	1.1	S	-	-
001529 BCNET16	1.6	S	1.1	H090100
001550 BCNET20	2	S	1.5	H090101
001530	2.4	S/M	1.8	H090208
001550A BCNET24	2.7	S/M	2.0	H090102
001551 BCNET30	3	S/M/L	2.3	H090203
001531	3.2	M/L	2.4	H090209
01551A BCNET35	3.6	M/L	2.7	H090104
001552 BCNET40	4	M/L	3.0	H090205
BCNET50	5	L	3.5	H090106S H090106L