

**ASIF TYPE SCREWS
HEXAGONAL HEAD**

Cortical CS Cortical Self Tap CSST Cancellous CAN	CORTICAL	CORTICAL	SELF-TAPPING CORTICAL	SMALL CANCELOUS	CORTICAL	SELF-TAPPING CORTICAL	SMALL CANCELOUS	CORTICAL	SMALL CANCELOUS	SELF-TAPPING CORTICAL	CANCELOUS	CORTICAL	CORTICAL	CANCELOUS
	1.5	2.0	2.0	2.2	2.7	2.7	3.0	3.5	3.5	3.5	4.0	4.5	5.5	6.5
JACOBS HARD PILOT DRILL SOFT	H090100	H090101 090101	H090101 090101	H090101 090101	H090102 090102	H090102 090102	H090102 090102	H090112 090112	H090102 090102	H090112 090112	H090102 090102	H090105S 090105	H090107 090107	H090105S 090105
	1.1	1.5	1.5	1.5	2.0	2.0	2.0	2.5	2.0	2.5	2.5	3.2	4.0	3.2
SYNTHES FIT	S090100S	S090101S	S090101S	S090101S	S090102S	S090102S	S090102S	S090112S	S090102S	S090112S	S090102S	S090105S	S090107	S090105S
TAP Quick Fit Tap 'T' Bar	TD0015 TT0015	TD0020 TT0020	N/A	N/A	TS00027 TD00027	N/A	N/A	TS12535 TT12535	TS17535 TT17535	N/A	TS17535 TT17535	TS00045 TT00045	TS00055	TS00065 TT00065
JACOBS HARD CLEARANCE DRILL SOFT	H090101 090101	H090102 090102	H090102 090102	H090102 090102	H090104 090104	H090104 090104	H090205	H090106S 090106	H090106S 090106	H090106S 090106	H090107 090107	H090108S 090108	H090109	6.5
	1.5	2.0	2.0	2.0	2.7	2.7	3.0	3.5	3.5	3.5	4.0	4.5	5.5	
SYNTHES FIT	S090101S	S090102S	S090102S	S090102S	S090104S	S090104S		S090106S	S090106S	S090106S	S090107	S090108S	S090109	
SCREW DRIVER	TYPE 1.5 HEX	1.5 HEX	1.5 HEX	1.5 HEX	2.5 HEX	2.5 HEX	2.5 HEX	2.5 HEX	2.5 HEX	2.5 HEX	2.5 HEX	3.5 HEX	3.5 HEX	3.5 HEX
	CODE SDD1500	SDD1500	SDD1500	SDD1500	SDH2540	SDH2540	SDH2540	SDH2540	SDH2540	SDH2540	SDH2540	SDH3500	SDH3500	SDH3500